

Study of Zeeman effect

Swaroop Ramakant Avarsekar*

School of Physical Sciences, National Institute of Science Education and Research, HBNI, Jatni -752050, India

(Dated: November 20, 2022)

I. INTRODUCTION AND THEORY

Zeeman effect is splitting of spectral lines under time independent magnetic field, due to interaction between the magnetic field and the spin angular momentum of the electrons where each of the electrons has its spin. When magnetic field is introduced, the degeneracy is broken and transition lines are now split so that there are more lines in emission spectra, which can be visualised by Fabry Perot's etalon resulting in circular fringes. The lines split in the pattern indicates Zeeman effect.




Figure 1: Splitting of components in magnetic field

If r is the radius of the lines, then

$$r_{n+1}^2 - r_n^2 = 2f^2/n_o = \text{constant} \quad (1)$$

This implies that

$$\Delta_a^{p+1,p} = r_{p+1,a}^2 - r_p^2, a = \Delta_b^{p+1,p} = r_{p+1,b}^2 - r_p^2, b \quad (2)$$

$$\therefore \Delta_a^{p+1,p} = \Delta_b^{p+1,p} \quad (3)$$

$$\delta_{a,b}^{p+1,p} = r_{p+1,a}^2 - r_{p+1,b}^2 \quad (4)$$

* swaroop.avarsekar@niser.ac.in

$$\Delta bar{\nu} = \delta / (2t\Delta) \quad (5)$$

where t is the distance between etalon plates.

The energy difference to calculate Bohr's magneton is as follows,

$$\Delta E = E_{L,M_L} - E_{L-1,M_{L-1}} = hc\Delta\bar{\nu}/2 \quad (6)$$

ΔE is proportional to the magnetic field as

$$\Delta E = \mu_B B \quad (7)$$

where μ_B is Bohr's magneton = $9.27 \times 10^{-4} J/T$

$$\mu_B = hc\bar{\nu}/2/B \quad (8)$$

μ_B can be obtained by linear fit between $\bar{\nu}$ and B with

II. EXPERIMENT

A. Apparatus




Figure 2: Schematic diagram for experimental setup

We would require Fabry-Perót etalon, permanent magnets, Cd lamp, optical bench, two lenses with $f=50$ mm and 300 mm, slide mounts, iris diaphragm to control amount of light, polarizing filter, CMOS camera with holder and a computer with respective software installed on it.

B. Procedure

Turn on the Cadmium light source which is placed between permanent pole pieces. The permanent magnets should be calibrated with variation of magnetic field by changing the distance between pole pieces. Align the etalon, lenses and camera such that a clear image is obtained from the computer screen. With varying the magnetic field the radius of lines varies. With the magnetic field turned on in the absence of the analyser three lines can be seen simultaneously in the normal Zeeman effect in transversal observation. In the case of the anomalous Zeeman effect three groups of three lines appear. Inserting the analyser in the normal Zeeman effect two σ -lines can be observed if the analyser is in the vertical position, while only the π -line appears if the analyser is turned into its horizontal position(transversal Zeeman effect). In the anomalous Zeeman effect there are two groups of three σ -lines in vertical polarization and one group of three π -lines in horizontal polarization. Turning the magnetic system by 90 degrees the light coming from the spectral lamp parallel to the direction of the field (longitudinal) can also be studied through the holes in the pole pieces. It can be shown that this light is circular polarized light (longitudinal Zeeman effect). A $\lambda/4$ -plate is generally used to convert linear into elliptical polarized light. In this experiment the plate is used in the opposite way. With the $\lambda/4$ -plate inserted before the analyser, the light of the longitudinal Zeeman effect is investigated. If the optical axis of the $\lambda/4$ -plate coincides with the vertical, it is observed that some rings disappear if the analyser is at an angle of +45 degrees with the vertical while other rings disappear for a position of -45 degrees . That means that the light of the longitudinal Zeeman effect is polarized in a circular (opposed way). The π -lines are longitudinally not observable.

C. Precautions

1. Always allow sufficient time for lamp to glow before taking readings
2. Do not place hands near the lamp source. It is very hot.
3. Do not touch optical surfaces with fingers.
4. Switch off the power supply before making any changes.

III. OBSERVATION AND ANALYSIS




Figure 4: Calibration curve for magnetic pole pieces.



Figure 3: Experimental setup in laboratory



Figure 5: Normal Zeeman effect



Figure 6: Longitudinal Zeeman effect




Figure 7: Anomalous Zeeman effect




Figure 8: Anomalous Zeeman effect in longitudinal direction with polariser.

Table I: Table for measurement radius of lines

Vernier Reading	Component	Radius			$r_1^2 (\mu\text{m}^2)$	$r_2^2 (\mu\text{m}^2)$	$r_3^2 (\mu\text{m}^2)$
		$r_1 (\mu\text{m})$	$r_2 (\mu\text{m})$	$r_3 (\mu\text{m})$			
41	σ_+	2998.88	4921.58	6226.76	8993281.254	24221949.7	38772540.1
	σ_-	1214.09	4155.9	5677.94	1474014.528	17271504.81	32239002.64
42	σ_+	2990.25	4971.72	6259.09	8941595.063	24717999.76	39176207.63
	σ_-	1517.98	4249.52	5715.25	2304263.28	18058420.23	32664082.56
44	σ_+	2983.79	4948.42	6187.39	8903002.764	24486860.5	38283795.01
	σ_-	1893.39	4362.2	5841.22	3584925.692	19028788.84	34119851.09
45	σ_+	2986.57	4917.36	6213.49	8919600.365	24180429.37	38607457.98
	σ_-	1996.95	4432.3	5862.09	3987809.303	19645283.29	34364099.17




Figure 9: Plot of wave number versus magnetic field

Table II: Calculation of Delta

Vernier Reading	$\delta (\mu\text{m}^2)$			δ
	δ^1	δ^2	δ^3	
41	7519266.726	6950444.886	6533537.454	7001083.022
42	6637331.782	6659579.528	6512125.066	6603012.125
44	5318077.072	5458071.656	4163943.924	4980030.884
45	4931791.062	4535146.08	4243358.812	4570098.651

Table III: Calculation of $\Delta\nu$ and Calibration of pole pieces.

Vernier Reading	B (T)	Component	$\Delta (\mu\text{m}^2)$		Δ	$\Delta\nu$ (m^{-1})
			$\Delta^{2,1}$	$\Delta^{3,2}$		
41	0.515	σ_+	15228668.44	14550590.4	15136061.74	50.0587921
		σ_-	15797490.28	14967497.83		
42	0.45	σ_+	15776404.7	14458207.87	15148607.96	47.21252557
		σ_-	15754156.95	14605662.33		
44	0.365	σ_+	15583857.73	13796934.52	14978929.41	35.60796663
		σ_-	15443863.15	15091062.25		
45	0.325	σ_+	15260829	14427028.61	15016036.87	32.67688978
		σ_-	15657473.99	14718815.88		

From this plot we obtain the slope as $98.26 \pm 13.31 (\text{mT})^{-1}$. With the slope obtained we use equation (8),

to obtain Bohr's magneton as $(9.76 \pm 1.32) \times 10^{-24} JT^{-1}$. The error can be calculated as

$$\frac{\delta\mu_B}{\mu_B} = \frac{\delta m}{m} \quad (9)$$

where m is slope, δ is associated error of that quantity.

The literature value of Bohr's magneton is $9.27 \times 10^{-24} JT^{-1}$ with the relative error is 4.97 %

IV. CONCLUSION AND SUMMARY

Zeeman effect is splitting of spectral lines under time independent magnetic field, due to interaction between the magnetic field and the spin angular momentum of the electrons where each of the electrons has its spin. Fabry Perot etalon was used to observed spectral lines. With the magnetic field turned on in the absence of the analyser three lines can be seen simultaneously in the normal Zeeman effect in transversal observation. In the case of the anomalous Zeeman effect three groups of three lines appear. The various Zeeman effects with

polariser was studied. We obtained Bohr Magneton's value $(9.76 \pm 1.32) \times 10^{-24} JT^{-1}$ with the relative error as 4.97 %.

The error can be reduced by properly adjusting the distance of lenses, camera and etalon such that the lines are visible sharp to facilitate the determination of radius. The apparatus should be in ambient temperature to avoid optical error artifacts. Random errors and human errors may also have contributed.

V. REFERENCES

1. https://www.niser.ac.in/sps/sites/default/files/basic_page/Zeemaneffect-manual.pdf
2. Zeeman effect old manual SPS NISER
3. https://en.wikipedia.org/wiki/Fabry%E2%80%93P%C3%A9rot_interferometer
4. https://www.nikhef.nl/~h73/kn1c/praktikum/phywe/LEP/Experim/5_1_10.pdf